

<b>CE Technical File</b> <b>First Aid Dressing Bandage</b>	File No: BT-CE-TF-1S-01_2.1
	Rev. No: A/2

# Instruction for use

Product Name: First Aid Dressing Bandage  
File No.: BT-CE-TF-1S-01\_2.1  
Version: A/1

Complied by:	<u>Echo Yan</u>	Date:	<u>2024.01.15</u>
Reviewed by:	<u>Beck Gao</u>	Date:	<u>2024.01.15</u>
Authorized by:	<u>Taylor Yang</u>	Date:	<u>2024.01.15</u>

**Bestreat Safety & First Aid Solution Co., Ltd.**

<b>CE Technical File</b>  <b>First Aid Dressing Bandage</b>	File No: BT-CE-TF-1S-01_2.1
	Rev. No: A/2

### Revision history

Revision Chapter	Revision Content	Revision Reason	Revision Date	Revised By	Document No. & Version	Remark
/	First issued	/	2022.09.15	/	BT-CE-TF-1S-01-2.1 A/1	/
Chapter II	Add the trade name\ BASIC UDI-DI/ Catalogue number	Incomplete information	2024.01.15	Echo Yan	BT-CE-TF-1S-01-2.1 A/2	/

#### 1.Product name

First Aid Dressing Bandage

Trade name: N/A

BASIC UDI-DI: 697570648100MT

Catalogue number: 200

#### 2. Brief introduction

The First Aid Dressing Bandage is a First Aid Dressing Bandage with aseptic bandage. Sterile First Aid Dressing Bandage material is polyesters or polyester and cotton blended, textile processing into rolls or tubular. The dressing is non-woven, and the wound surface is covered with a layer of porous material with non-absorptive properties and can penetrate the exudate. The dressing contacts the wound directly.

#### 3. Intended Use

It is intended for non-chronic wound care, or for wound dressings or limbs to provide binding force, to play a binding, fixed role.

#### 4. Models of product

##### 4.1 Model

According to the dressing material of Sterile First Aid Dressing Bandage, it can be divided into one type of non-woven fabric. Elastic bandages, Non-elastic bandages, With or without absorbent pad; Round shaped, Triangular shaped, Square shaped, Trauma bandages, etc.

##### 4.2 Specification:

<b>Specifications</b>	<b>Dressing</b>	<b>Bandage</b>	<b>The distance between the dressing and</b>
-----------------------	-----------------	----------------	--

<b>CE Technical File</b>  <b>First Aid Dressing Bandage</b>	File No: BT-CE-TF-1S-01_2.1
	Rev. No: A/2

					the shorter end of the bandage
	<b>b2</b>	<b>l2</b>	<b>b1</b>	<b>l1</b>	<b>l3</b>
	<b>±0.5</b>	<b>/</b>	<b>±0.5</b>	<b>Elastic±2 0</b>	<b>±3</b>
5cm x 225cm	5	5~8	5	225	15
5cm x 300cm	5	5~8	5	300	15
5cm x 400cm	5	5~8	5	400	15
5cm x 450cm	5	5~8	5	450	15
6cm x 225cm	6	6~10	6	225	15
6cm x 300cm	6	6~10	6	300	15
6cm x 400cm	6	6~10	6	400	15
6cm x 450cm	6	6~10	6	450	15
7.5cm x 225cm	7.5	7.5~10	7.5	225	15
7.5cm x 300cm	7.5	7.5~10	7.5	300	15
7.5cm x 400cm	7.5	7.5~10	7.5	400	15
7.5cm x 450cm	7.5	7.5~10	7.5	450	15
8cm x 225cm	8	8~10	8	225	15
8cm x 300cm	8	8~10	8	300	15
8cm x 400cm	8	8~10	8	400	15
8cm x 450cm	8	8~10	8	450	15
10cm x 225cm	10	10~15	10	225	15
10cm x 300cm	10	10~15	10	300	15
10cm x 400cm	10	10~15	10	400	15
10cm x 450cm	10	10~15	10	450	15
15cm x 225cm	15	10~20	10	225	15
15cm x 300cm	15	10~20	10	300	15
15cm x 400cm	15	10~20	10	400	15
15cm x 450cm	15	10~20	10	450	15

## 5.Contraindications

Nonwoven allergies are forbidden to use

## 6.Sterilization

Steam sterilization or Ethylene oxide sterilization (EO sterilization)

## 7.Precautions

- 1) This product is a disposable product.
- 2) Do not use expired products
- 3) Damaged inner packing is strictly prohibited.
- 4) When dressing, the blood circulation of body tissues should not be blocked due to excessive tightness.
- 5) When the bandage is used, it should be flattened to avoid reducing the bandage width.
- 6) The used bandages should be put into the medical dustbin.

<b>CE Technical File</b>  <b>First Aid Dressing Bandage</b>	File No: BT-CE-TF-1S-01_2.1
	Rev. No: A/2

## 8. Instruction for use

- 1) Tear off along the upper edge of the inner package to remove the product.
- 2) Open the bandage at one end of the bandage until the patch is exposed.
- 3) The center of the dressing is aligned with the wound surface, and after the end of the bandage is fixed, a circular continuous wrapping bandage is used.
- 4) Fix the other end of the bandage with a medical tape.
- 5) Wounds can be treated with iodophor, alcohol, etc., and bandage is used.

## 9. Major raw materials

Polyester fiber and dressing (The dressing is called a water absorbent pad. Its main raw material is nonwovens)

The chemical composition:

Bandage: Polyester fiber 100%

Absorbent pad: Polythene 15% + Viscose 60% + Polyester fiber 25%

## 10. Storage conditions

Should be stored in a well-ventilated room without corrosive gas.







## 11. Transport Conditions

Sunscreen and rainproof during transportation








## 12. shelf life of product

5 years and 3 months

## 【Definitions of graphics, symbols, abbreviations】

	DO NOT REUSE		USE BY
	BATCH CODE		" CE " Marking
	DATE OF MANUFACTURE		AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY

<b>CE Technical File</b>  <b>First Aid Dressing Bandage</b>	File No: BT-CE-TF-1S-01_2.1
	Rev. No: A/2

	MANUFACTURER	<b>MD</b>	Medical devices
	KEEP AWAY FROM SUNLIGHT		KEEP DRY
	STERILIZED USING ETHYLENE OXIDE		STERILIZED USING STEAM OR DRY HEAT
	DO NOT USE IF PACKAGE IS DAMAGED		SINGLE STERILE BARRIER SYSTEM

### 【Manufacturer Information】



#### Manufacturer

**Company Name:** Bestreat Safety & First Aid Solution Co., Ltd

Address: Room 401, 4th Floor, Building 2, Bonded Zone, No.378 Gexin Avenue, Dongxihu District, Wuhan, 430040 Hubei, P.R. China

Contact person: Echo Yan

E-mail: qm@bestreatkit.com

### 【Authorised Representative Information】



#### Authorised representative in EEC

Name: MedPath GmbH

Add: Mies-van-der-Rohe-Strasse 8,80807 Munich, Germany

Tel: +49(0)89 189174474

Dimdi Code: DE/0000047823

E-mail: info@medpath.pro