



Benefits >>>

Effectively reduce the lateral tension on both sides of the wound, prevent the scar from widening and hyperplasia at the suture needle; effectively dilute the wound and reduce the chance of scar formation.

Application >>>

It is used to provide binding force on wound dressings or limbs, so as to bandage and fix.

Product Structure >>>

The product is composed of zipper components, acrylic vinegar adhesive, polyurethane film, release paper, etc.; when in use, the polyurethane film is fixed on the skin surface on both sides of the wound, and when the zipper component is pulled, the skin is pulled toward the center of the wound. Move closer to close the wound, and then through the locking force of the spines, the wound is in a closed state to achieve the purpose of suturing.

PRODUCT ADVANTAGE

- Through physical action, it reduces the tension around the wound, and has a unique design of dynamic compression and decomposition of tension, which can reduce the risk of scar formation.
- It can be used for the combination of common epidermal wounds in daily life to avoid the secondary injury of suture needles.
- Individually packaged – Different from other clean packaged products, this product is treated with ethylene oxide and can be directly pasted on the wound to avoid cross-infection.
- Special adhesive, firmly fixed, waterproof and breathable, low skin irritation, low allergy;
- The hollowed-out design in the center can observe the wound healing at any time without affecting dressing changes after surgery.
- It can meet different radian wounds and is widely used in many types of wounds.