

# YJ-UR200

## Instructions

1. Test species: humans and animals
2. Test categories: water and urine
3. Test items:  
(Urine test 14 parameters, water test 14 parameters)
4. Test principle: super-high lumimosity cold light
5. Screen type: 7-inch capacitive touch screen.
6. Test speed: up to 120 tests/hour at continuous test mode
7. Use test paper: urs-11-iii, h12ma, w-14, urs-10t, urs-14
8. Data storage: more than 1000000 results:
9. Printing mode: thermal printer
10. External interface: microsd (tf), maximum support 8g, usb2.0 (2)
11. Data transmission: standard rs232c (115200bps)
12. Network connection: 10/100/1000mbps ethernet, wifi
13. Input device: keyboard and mouse
14. Use language: support chinese, english, russian and many other national languages.
15. Power supply: 100-240v adaptive ac power adapter
16. Net weight: 3kg
17. Volume: 355\*300\*145mm



## Product Features

1. Adopt with linux operation system and potable human-machine interaction interface.
2. High speed cpu, main frequency 528mhz
3. Online upgrade system
4. Large capacity data storage
5. External serial /usb bar code reader
6. Wireless connection network
7. Built in thermal printer (remove "equipped with an external pln printer")
8. Program memory 64m+512k, running memory 32m+96k  
The single chip control of the whole system enhances the stability and reliability of the system.
9. Built in thermal printer can print two sizes, and the definition is adjustable.



Henan Forever Medical Co.,Ltd.

24

in-vitro diagnostic medical device (IVD)

## Instructions

- \* Test Principle: super-high lumimosity cold light source reflection
- \* Test parameters: VC, LEU, NIT, URO, PRO, PH, BLD, SG, KET, BIL, GLU
- \* Test speed: Up to 120 tests/hour at Continuous Test Mode
- \* Memory Size: 1000 test results ; 8 , 10, 11 compatible (optional)
- \* Absord the excess urine automatically, prevent from cross-inflection of urine
- \* Thermal Printer Inside, with outer stylus Printer interface
- \* PC connected with RS232
- \* Data Communication Interface Baud Rate: 9600bps
- \* Language: English and Chinese
- \* Environment Condition: 0-40°C, RH < 85%



# YJ-UR200B

## Urine Analyzer

HENAN FOREVER MEDICAL CO.,LTD.  
Room 806, Floor 8, Wanda Center, Jinshui District, Zhengzhou, Henan Province, China  
FAX: 86-371-61315851  
Email: judy@yjcompany.com  
Whatsapp/TEL.: 0086-15290808151  
Web: www.forermedchina.com



Alibaba Shop



Made-In-China Shop



WEBSITE