

CareFix®

HAND , FOOT & ANKLE

Internal Fixation System

Comprehensive technical Solution For Trauma and Orthopedics



Shanghai CareFix Medical

Cannulated Screw System

Cannulated screws have been widely used in hand and foot & ankle surgery, especially for small bone fractures, joint fusion, and osteotomy.

Headless Compression Cannulated Screw

- Φ2.5x10-30/2mm
- Φ3.0x12-40/2mm
- Φ3.5x14-40/2mm
- Φ4.0x16-50/2mm
- Φ4.5x20-50/2mm, 55-60/5mm
- Φ6.5x40-50/2mm, 55-120/5mm



Double Threaded Compression Cannulated Screw

- Φ2.0x10-20/2mm
- Φ2.5x10-30/2mm
- Φ3.0x10-40/2mm
- Φ4.0x16-50/2mm, 55-100/5mm
- Φ4.5x20-50/2mm, 55-70/5mm
- Φ6.5x40-120/5mm (Long/Short Thread)




Compression Cannulated Screw

- Φ3.5x14-50/2mm
- Φ4.5x30-50/2mm, 55-70/5mm
- Φ6.5x40-120/5mm(16,32mm Thread)



Full Threaded Cannulated Screw

- Φ3.5x18-50/2mm, 60-100/10mm
- Φ4.5x20-50/2mm, 50-60/5mm, 70-110/10mm




Washer

- Φ3.2x1.2mm
- Φ4.3x1.5mm
- Φ6.3x1.5mm



Subtalar joint stabilizer

Subtalar joint stabilizer I



Subtalar joint stabilizer II




Subtalar joint stabilizer III



Advantages:


1. Bionic design: Reduced pain, low fall-off
2. Comprehensive range: Covers major standards
3. Minimally invasive: Rapid recovery
4. Ergonomic tools: Simplified operation

Metallic Bone Pins



1. Main Categories: Kirschner wire, double-headed kirschner wire, olive pin, steinmann pin, snap-off screw
2. Specification: 0.8, 1.0, 1.2, 1.5, 1.8, 2.0, 2.5, 3.0, etc
3. Material: Stainless steel and titanium alloy
4. Sterilization and non-sterilized packaging

Snap-off Screw



U-Screw

Round U-Screw

10x15x1.5
13x20x1.5

Anti-Slip U-Screw

10x9x1.5
13x11x1.5
15x13x1.5

Compression U-Screw

15x13x3.0
18x15x3.0
20x15x3.0



Indications:

- ▶ Fixation of metaphyseal and intra-articular fractures
- ▶ Internal fixation of pelvic, flat, irregular, and small bones
- ▶ Auxiliary fixation of long bone fractures
- ▶ Hand or foot fracture, osteotomy and joint fusions
- ▶ Epiphyseal block in children

Material: Stainless Steel

Micro Compression Locking Plate 1.5/2.0

Limited Compression Locking Plate 4-8 holes



LCP Reconstruction Plate 10 holes



LCP Reconstruction (head 2 holes) 3-5 holes



LCP Reconstruction (head 3 holes) 3-5 holes



LCP Reconstruction (head 2 holes) 3-5 holes



LCP Reconstruction (head 3 holes) 3-5 holes



Mini Compression Locking Plate 2.4/2.7

Limited-Compression Locking Plate 4-8 holes



Reconstruction Compression Locking Plate 12 holes



Reconstruction Compression Locking Plate (head 2holes) 3-7 holes



Reconstruction Compression Locking Plate (head 3holes) 3-7 holes



Reconstruction Compression Locking Plate (head 2holes) 3-7 holes




Reconstruction Compression Locking Plate (head 2holes) 3-7 holes




Features:

Reconstruction design: easily malleable or shearable. Anatomical design: enhanced bone fit.


Limited contact: minimizes periosteal blood flow disruption. Low notch design: minimizes soft tissue irritation. LCP combi hole: optional compression or angular locking.



1.5Cortical Screw 6-20mm/2mm
2.0Cortical Screw 8-24mm/2mm



1.5Locking Screw 6-20mm/2mm
2.0Locking Screw 8-24mm/2mm



2.4Cortical Screw 8-30mm/2mm
2.7Cortical Screw 10-30mm/2mm




2.4Locking Screw 10-30mm/2mm
2.7Locking Screw 10-40mm/2mm

Hind Foot and Ankle Locking Plate System

Minimally Invasive Calcaneal Locking Plate 3.5/4.0 Left/Right




Calcaneal Locking Plate III 3.5/4.0 (Small,Medium,Large) Left/Right



LCP Calcaneus Plate II 3.5/4.0 (Small,Medium,Large) Left/Right



LCP Calcaneus Plate V 3.5/4.0 (Small,Medium,Large) Left/Right



LCP Distal Fibular Plate 3.5/4.0 3-10Holes



LCP Medical Distal Tibia Plate 3.5/4.0
5/7/9/11/13 holes Left/Right




LCP Lateral Distal Tibia Plate 3.5/4.0
5/7/9/11/13 holes Left/Right



Anterior Dista Tibial Locking Plate 3.5/4.0 3/5/7/9 holes





Foot and Ankle Special Screws












3.5mm low-profile locking screw
12-50mm/2mm



3.5mm low-profile cortical screw
12-40mm/2mm




Foot and Ankle Orthopedic Locking Plate

| Picture | Product Name | Specification | Length |
|---|-------------------------------|---------------|----------------------|
|  | 1st MTP Fusion LCP 2.4/2.7 | 5° | 52mm Left/Right |
| | | 10° | 57mm Left/Right |
|  | T-type Fusion LCP 2.4/2.7 | 2 holes | 35mm |
| | | 2 holes | 42mm |
| | | 4 holes | 60mm |
|  | L-type Fusion LCP 2.4/2.7 | 2 holes | 37mm Left/Right |
| | | 2 holes | 44mm Left/Right |
| | | 4 holes | 62mm Left/Right |
|  | Cloverleaf Fusion LCP 2.4/2.7 | 2 holes | 38mm |
| | | 2 holes | 45mm |
| | | 4 holes | 64mm |
|  | Locking plate 2.4/2.7 | 4 holes | 22mm |
| | | 4 holes | 24mm |
| | | 4 holes | 30mm |
| | | 4 holes | 36mm |
|  | Step Locking plate 2.4/2.7 | 4 holes | Step 0,1,2,3,4,5,6mm |
| | Step Locking plate 3.5/4.0 | | Step 2,4,6,8,10mm |
|  | Opening Wedge 2.4/2.7 | 4 holes | Opening Wedge 2-7mm |
| | Opening Wedge 3.5/4.0 | | Opening Wedge 6-10mm |
|  | W-type 2.4/2.7 | 4-10 holes | 21-57mm |
|  | Navicular | 11 holes | / |

Ankle Fusion Locking Plate

Multi-Functional Ankle Fusion Plate




China's First: 6-Piece ankle fusion locking plate.

Anatomical design, cross-joint fixation

Multiple Approaches: Anterior, Lateral, Posterior

Specialized compression & Spreader Forceps

Indications: Arthritis, Trauma (12 years old or above)




Anterior Ankle Fusion Locking Plate



Anterior Ankle Fusion Locking Plate



Posterior Ankle Fusion Locking Plate



Lateral Ankle Fusion Locking Plate

Distal Tibia Osteotomy plate

Distal Tibia Osteotomy T-Plate 3.5/4.0
3/5/7 holes



Distal Tibia Osteotomy Locking L-Plate 3.5/4.0
4/5/6 holes Left/Right



Ankle Fusion Intramedullary Nail

Intramedullary Nail (IM nail) :



| Diameter | About | Length |
|----------|------------|-----------------|
| 9mm | left/right | 150、 200、 250mm |
| 10mm | left/right | 150、 200、 250mm |
| 11mm | left/right | 150、 200、 250mm |
| 12mm | left/right | 150、 200、 250mm |


Interlocking Screw :



| Diameter | Length |
|----------|-----------------------|
| 5.0 | 25-90 (5mm interval) |
| 6.0 | 25-120 (5mm interval) |




Foot & Ankle Basic Surgical Instruments



Reamer


---Specialized Design and Production



Distraction Piler



Foot & Ankle Basic Instrument Kit



Foot & Ankle Locking Plate Instrument Kit



Micro Locking Plate System Instrument Set

Shanghai CareFix Medical Instrument Co.,Ltd

Web: www.carefix-ortho.com Tel: +86 400 098 2859

Email: carefixsh@gmail.com WhatsApp: 18918690107

Add: No.31 Lane 1515, Rongle East Road, Songjiang District, Shanghai, China

