



## YJ-FA516 FA526 SYRINGE PUMP



### Features:

- 3.6-inch LCD screen, better visualization
- Wide range of flow rate: 0.1-1500ml/h
- Numeric keypad design, easy-to-set
- Automatic recognition of syringe specifications (10-60ml)
- Comprehensive alarm function (360 visible & audible)
- Multiple anti-bolus detection technologies
- Integrated design (FA526), better bedside space efficiency
- Optional: WiFi module





Technical Specifications				
Model	YJ-FA516		YJ-FA526	
Channel	Single		Double	
Flow Rate Range(ml/h)	10ml:0.1-300ml/h; 20ml:0.1-600ml/h 30ml:0.1-900ml/h 50/60ml:0.1-1500ml/h Step by 0.1			
Accuracy	±2%			
Volume Range(ml)	0-9999.9 Step by 0.1			
Purge/Bolus Rate	10ml:0.1-300ml/h; 20ml:0.1-600ml/h 30ml:0.1-900ml/h 50/60ml:0.1-1500ml/h Step by 0.1			
Blocking alarm value	High	20ml,30ml,50ml:90Kpa±30Kpa 10ml:100Kpa±30Kpa		
	Low	20ml,30ml,50ml:50Kpa±30Kpa 10ml:70Kpa±30Kpa		
KVO Rate	0.1ml/h ±0.02 ml/h			
Purge / Bolus Function	√			
Alarm	Occlusion, Pause Overtime, Low Battery, Near End, VTBI Completion, Syringe Near Empty, Syringe Empty, Incorrect Syringe Specification, Syringe Dislocation, Motor Malfunction, Battery Exhausted, Battery Malfunction.			
Pause Overtime Alarm	1min50s-2min			
Power & Battery	~100V~240V 50Hz/60Hz DC:12V Inner Ni-MH battery 12V			
Weight	2.03 Kg		3.27 Kg	
Battery Time	Ni-MH battery support 7 hours for single channel, 5.5 hours for double channel at intermediate rate under full of charge.			
Classification	Class I, type CF,IPX4			
Size	290*105*138mm		360*195*138mm	